INSTRUCTIONS FOR USE

For use outside USA only.

INTEGRA® DERMAL REGENERATION TEMPLATE

Dermal Regeneration Template.
Bilayer membrane system used as a template
for dermal regeneration.

INTEGRA® DERMAL REGENERATION TEMPLATE SINGLE LAYER

Dermal Regeneration Template Single Layer.
Single layer membrane system used as a
template for dermal regeneration.

English

STERILE R

DESCRIPTION

INTEGRA® Dermal Regeneration Template is a bilayer membrane system for skin replacement. The dermal replacement layer is made of a porous matrix of fibers of cross-linked bovine tendon collagen and a glycosaminoglycan (chondroitin-6-sulfate) that is manufactured with a controlled porosity and defined degradation rate. The epidermal substitute layer is made of a thin polysiloxane (silicone) layer to control moisture loss from the wound.

INTEGRA Dermal Regeneration Template Single Layer consists only of the dermal regeneration layer and is available to add extra thickness to the dermal regeneration layer, when deep wounds are to be treated.

INTEGRA template is provided sterile. The inner foil pouch and product should be handled using sterile technique. INTEGRA template should not be re-sterilized.

INDICATIONS

INTEGRA Dermal Regeneration Template is indicated for the postexcisional treatment of full-thickness and partial-thickness injuries where sufficient autograft is not available at the time of excision or not desirable due to the physiological condition of the patient. INTEGRA Dermal Regeneration Template is also indicated for use in reconstruction of postexcisional, full-thickness defects of the integument where there is, in the opinion of the treating surgeon, a potential benefit to the patient by improving the reconstructive outcome or decreasing their mortality/morbidity.

CONTRAINDICATIONS

Use of INTEGRA Dermal Regeneration Template is contraindicated in patients with known hypersensitivity to bovine collagen, chondroitin or silicone materials.

INTEGRA template should not be used on clinically diagnosed infected wounds.

SAFETY

INTEGRA Dermal Regeneration Template is manufactured from collagen obtained from bovine deep flexor tendon, which is classified by European Standards as Category C material (no detectable infectivity for Bovine Spongiform Encephalopathy (BSE)). Bovine Tendon is known to be one of the purest sources of Type I collagen that is commercially available.

The collagen used to manufacture INTEGRA Dermal Regeneration Template is currently used in the manufacture of artificial skin, absorbable hemostatic sponges, and absorbable wound dressings. The manufacturing process for INTEGRA Dermal Regeneration Template meets USA and European Standards for animal tissue sourcing, handling and inactivation of viruses and transmissible agents.

This process involves a treatment with sodium hydroxide that is a recognized method of inactivation of Spongiform Encephalopathy pathogens.

A viral inactivation study for the INTEGRA Dermal Regeneration Template manufacturing process was conducted by an independent certified laboratory. In this study, the sodium hydroxide reduced the viral titer to non-detectable levels for the following viral strains: Human Immunodeficiency Virus Type 1 (HIV), Bovine Viral Diarrhea (BVD), Infectious Bovine Rhinotracheitis (IBR), Parainfluenza Virus Type 3 (PI3), Vesicular Stomatitis (VSV).

WARNINGS

Excision of the wound must be performed thoroughly to remove all coagulation eschar or nonviable tissue from the wound bed. INTEGRA Dermal Regeneration Template will not take to nonviable tissue. Leaving any remaining nonviable tissue may create an environment for bacterial growth.

Hemostasis must be achieved prior to applying INTEGRA template. Inadequate control of bleeding will interfere with the incorporation of INTEGRA template.

PRECAUTIONS

There have been no clinical studies evaluating INTEGRA Dermal Regeneration Template in pregnant women. Caution should be exercised before using INTEGRA template in pregnant women. Such use should occur only when the anticipated benefit clearly outweighs the risk.

In clinical trials, the use of INTEGRA template was evaluated in a small number of patients with chemical, radiation, or electrical burns. A surgeon's decision to use INTEGRA template on these wounds should be based on their evaluation of the wound and its suitability to excisional therapy, the likelihood that a viable wound bed will be created by excision, and whether the possible benefit outweighs the risk in this patient population.

INTEGRA template should be applied on the day of excision. Delaying the application of INTEGRA template may substantially impair the take of the material.

Appropriate techniques to minimize pressure and shearing should be used to reduce risk of mechanical dislodgement, especially if the product is applied over joints or on pressure areas such as the back or buttocks.

Placing the patient in hydrotherapy immersion may interfere with proper incorporation of INTEGRA Dermal Regeneration Template and cause premature separation of the silicone layer and nonadherence of INTEGRA template.

Caution must be employed not to remove the newly formed neodermal tissue when removing the silicone layer. The neodermis must NOT be excised off the wound.

ADVERSE EVENTS

INTEGRA Dermal Regeneration Template has been found to be well tolerated in 4 prospective clinical trials involving 444 burn patients, 1,200 wound sites and reconstructive surgery studies involving 119 patients and 157 wound sites. There were no reports of clinically significant immunological or histological responses to the implantation of INTEGRA template. There were no reports of rejection of INTEGRA template in the clinical evaluations.

There were adverse events reported in the INTEGRA template clinical trials which are those commonly reported in patients with severe burns. Adverse events reported in the burn and reconstructive surgery clinical evaluations included infection, failure to incorporate into the wound bed, partial or complete graft loss, hematoma, fluid under the silicone layer, shearing or mechanical loss of the product, separation of the silicone layer and delayed healing.

In the clinical trials, data were collected regarding wound colonization or infection. The consequences of infection at sites treated with INTEGRA template included partial or complete loss of take of INTEGRA template. Infection rates in burn wound sites treated with INTEGRA template in the clinical trials ranged from 14 to 55%. In the Post approval study involving 216 patients and 841 wound sites, overall infection rate was 16.3%.

SINGLE-USE DEVICE

INTEGRA Dermal Regeneration Template is supplied in a single-use package and is guaranteed to be sterile and non-pyrogenic unless opened or damaged. The product is intended for use as an absorbable implant and is not to be reused. Reuse of the device can result in contamination and/or disease transmission. Any attempt to resterilize or reuse the product/components will damage the matrix and impair its ability to function as intended. All unused pieces must be discarded.

INFORMATION FOR USE

INTEGRA Dermal Regeneration Template facilitates the formation of a neodermis by the body. The collagen/glycosaminoglycan dermal portion serves as a template for the infiltration of fibroblasts, macrophages, lymphocytes, and capillary endothelial cells, which form the neovascular network. As healing progresses, the collagen-GAG layer is resorbed and new collagen is deposited by the fibroblasts to form the neodermis. Upon adequate vascularization of the neodermis and the availability of the donor autograft, the silicone layer is removed and a thin, meshed or unmeshed layer of epidermal autograft is placed over the neodermis. Cells from the epidermal autograft grow and form a mature epidermis thereby closing the wound, and resulting in a functional dermis and epidermis.

Histological Evaluation of INTEGRA Dermal Regeneration Template

Three hundred thirty-six serial biopsies were obtained from 131 patients participating in the multicenter clinical trial ranging from 7 days to 2 years after application of INTEGRA template. A histological study of the wound healing in the burned areas was conducted. An intact dermis was achieved with regrowth of apparently normal reticular and papillary dermis. No scar formation appeared in the biopsies of patients examined.

Patient Evaluation and Selection of Sites

In patients receiving INTEGRA Dermal Regeneration Template for injuries, the extent of the defect that requires prompt excision and grafting is identified.

Patients undergoing reconstructive procedures should have appropriate preoperative evaluations.

INTEGRA template may be applied to all excised wound sites.

SURGICAL APPLICATION OF INTEGRA DERMAL REGENERATION TEMPLATE

Scheduling Surgery

INTEGRA template must be applied to a viable wound bed following surgical excision. For acute cases surgery may be scheduled as soon as the patient is stabilized, generally within 1 to 3 days after injury. Surgery may be staged as appropriate.

In reconstructive procedures surgery is scheduled according to the judgement of the surgeon.

Perioperative Antibiotics

It is recommended that perioperative antibiotics be administered according to the clinical judgment of the practitioner.

Product Preparation

INTEGRA Dermal Regeneration Template

1. Always handle INTEGRA template using aseptic technique.
2. Peel open the outer pouch and remove the inner foil pouch using sterile technique.
3. Place foil pouch flat on a sterile surface and peel it open.
4. Remove product, including the protective polyethylene cover sheets.
5. While holding the product with the tab remove one polyethylene cover sheet. Turn the product and remove the second polyethylene cover sheet.
6. Using the tab, the product can now be placed into a basin containing sterile saline solution. Carefully remove the tab from the product.
7. Rinse the product by immersion in sterile saline for 1-2 minutes.
8. Keep product in the basin until application.

INTEGRA Dermal Regeneration Template Single Layer

1. Always handle INTEGRA template using aseptic technique.
2. Peel open the outer pouch and remove the inner foil pouch using sterile technique.
3. Place foil pouch flat on a sterile surface and peel it open.
4. Remove product including the protective polyethylene sheets.
5. Separate product from the polyethylene sheets by slowly and carefully peeling back from one of the corner edges. Handle with care as product is fragile.
6. Place product into basin containing sterile saline solution.
7. Rinse the product by immersion in sterile saline for 1-2 minutes.
8. Keep product in the basin until application.

Meshing INTEGRA

INTEGRA template can be meshed at 1:1 ratio before the application, but must not be expanded. Meshing may improve the ability of INTEGRA template to conform to irregular surfaces and may improve take on exuding wounds.

The use of topical antimicrobials is recommended if INTEGRA template is meshed.

Wound Excision

Excision must be made to the level of viable tissue and meticulous hemostasis must be achieved before application of INTEGRA Dermal Regeneration Template.

Excisional techniques for INTEGRA template sites can be fascial, sequential, or tangential. It is absolutely critical to the successful take of INTEGRA template that excision is complete and that no devitalized tissue remains. It is recommended that cultures be taken (swab or quantitative) to monitor all areas.

When excising scar or releasing contractures, care should be taken to achieve a complete release including deep contractures.

Complete hemostasis must be achieved before application of INTEGRA template. The presence of hematoma will cause loss of INTEGRA template in the affected area. Broad area cauterization that could decrease wound bed viability should be avoided.

Shaping INTEGRA Dermal Regeneration Template to Fit the Wound

INTEGRA template should be shaped accurately to fit the excised wound margins to minimize scarring at these margins. INTEGRA template should not be overlapped onto unexcised burn eschar. It is easily cut with sterile scissors by placing the sheet of INTEGRA template over the open area and cutting exactly to the edge.

Applying INTEGRA to the Wound

It is critical that the collagen template layer be in direct contact with the excised wound. The silicone layer (identified by the black threads) must be placed out (away from the wound bed). Do not apply upside down; the black threads must be clearly visible.

INTEGRA Dermal Regeneration Template Single Layer should be placed in the wound bed in a similar fashion. INTEGRA Dermal Regeneration Template is then placed on top as the outer layer. Ensure that there are no folds and gaps between INTEGRA template and the wound bed or between the layers of INTEGRA template.

For optimal cosmetic results, place the INTEGRA template so that the suture lines between INTEGRA template sheets lie in Langer's lines. This will minimize the final appearance of the suture lines.

The INTEGRA template sheets are secured by staples or sutures placed in an interrupted fashion. If the INTEGRA template has been meshed, care must be taken not to spread the mesh. Care should be taken to achieve a primary closure between INTEGRA template and adjacent unburned skin or between sheets of INTEGRA template. Suture or staple each sheet of INTEGRA template in place independently. The INTEGRA template sheets may also be affixed to one another to assure that there is no gap. Adjust the area to ensure that there is no undue tension on an individual piece of INTEGRA template. The material should readily adhere and conform to the wound surface. Any air bubbles should be carefully removed by moving them to the edge of the sheet. Dressings should be applied over INTEGRA template according to the protocol of the practitioner.

Postoperative Care

Postoperative care, like that used following treatment with full sheet or meshed autograft, should include monitoring for hematomas, wound infection and patient sepsis.

The outer dressings should be changed as necessary. However, the inner dressing need not be disturbed unless there are problems requiring intervention. The attachment of the silicone layer should be examined. An antibacterial dressing may be used or the outer dressing can be soaked in an appropriate antimicrobial solution.

There must be no hydrotherapy immersion of the patient following INTEGRA template application while the silicone layer is in place.

Mechanical dislodgement of INTEGRA template should be avoided. Ambulation and physical therapy can be instituted according to the condition of the patient and judgement of the practitioner. All INTEGRA template sites must be securely covered with dressings before ambulation and/or physical therapy.

Staples or sutures should remain in place until the time of epidermal autografting. The staples or sutures help secure the silicone layer and decrease the likelihood of premature silicone layer separation.

EPIDERMAL AUTOGRAFTING

Identifying the Neodermis

The area of INTEGRA Template take (neodermis formation) should be evaluated before application of the epidermal autograft. Neodermis may be recognized by a yellow-orange color with occasional areas of light red. The neodermis should be firmly attached to the underlying tissue. The silicone layer should easily separate from the underlying neodermis.

Removal of Silicone Layer for Epidermal Autograft

The silicone layer of INTEGRA template may be removed when the collagen layer has been replaced by neodermis, usually 14 to 28 days after application of INTEGRA template. Time for neodermis formation may depend on the patient, anatomic location and cause of the defect.

Removal of the silicone layer and grafting may take place immediately after formation of the neodermis, if epidermal autograft is available. The removal of the silicone layer may be postponed until donor sites for epidermal autograft area are available.

The clinician must be careful when removing the silicone layer. The silicone layer can usually be removed using only forceps. Generally, it should peel off easily. Difficulties in removal may indicate that neodermis formation is incomplete. However, if the silicone is difficult to remove, a forceps and scalpel may be used to gently separate the silicone layer from the neodermis. Caution must be employed to not remove the newly formed neodermis tissue when removing the silicone layer.

Harvesting and Preparation of the Epidermal Autograft

Epidermal autograft can be taken from sites unsuitable for conventional autograft, for example small areas from which a large intact sheet would be impossible. If possible, the area should be matched for color and type of skin. The epidermal autograft should be at a thickness just sufficient to provide punctate bleeding of the donor site, approximately 0.15 mm (0.006 inches) (range 0.10mm-0.20 mm / 0.004-0.008 inches). Dermal tissue is not needed in the epidermal autograft, and should be minimized.

The thin epidermal autograft can be meshed up to a 4:1 ratio. The meshed epidermal autograft may be fragile and care should be taken in handling the graft.

Application of the Epidermal Autograft to the Neodermis

The epidermal autograft should be placed over the neodermis by spreading the meshed autograft. It should be spread as evenly as possible over the neodermis without leaving large open areas.

Completion of the epidermal autograft procedure should follow the standard protocol for full sheet or meshed autograft. The epidermal autograft should be anchored by sutures or staples. The dressing over the epidermal autograft should be similar to that used over conventional meshed autograft.

Postoperative Epidermal Autograft Care

Postoperative care of the thin epidermal autograft should follow the protocols for conventional sheet or meshed autograft.

Treatment as an Outpatient

INTEGRA Dermal Regeneration Template may be placed on patients who will be treated as outpatients after surgery. The above protocol should be followed and an appropriate wound bed maintained. The patient should be under the supervision of the treating physician. Care must be taken not to dislodge the INTEGRA template prior to epidermal grafting.

Epidermal grafts are initially fragile and the wounds should have protective dressings in place. This can generally be accomplished with wound maintenance and education of patient, family, home health nursing, and outpatient clinic personnel.

POTENTIAL POSTOPERATIVE PROBLEMS

Wound Colonization or Infection

Wounds having excessive discharge may require more frequent dressing changes and may require the use of appropriate antimicrobial intervention. After a diagnosed infection is controlled, new INTEGRA template or, if a part of the neodermis remains, thin epidermal autograft may be applied to the affected area.

Patient Sepsis

The dressings should be removed and wound sites (INTEGRA template or autograft) inspected for infection. Appropriate diagnostic and therapeutic procedures should be followed.

Hematoma

Areas of hematoma should be monitored and aspirated or excised as required. New INTEGRA template or autograft may be applied to the excised sites after bleeding is controlled.

Separation of INTEGRA Dermal Regeneration Template From the Wound Bed or Poor INTEGRA Dermal Regeneration Template Take

If INTEGRA template becomes separated from the wound bed prior to neodermis formation, remove the silicone layer and the attached dermal matrix and carefully inspect the wound bed for residual nonviable tissue or infection. If there is no nonviable tissue or infection, the wound may be regrafted with INTEGRA template after superficial debridement.

Areas of poor INTEGRA template take may be treated by reapplication of INTEGRA template or by application of conventional autograft.

Fluid Accumulation and Premature Silicone Layer Separation

Fluid accumulation or premature silicone layer separation must be treated to prevent infection or granulation tissue. Small areas of fluid accumulation under the silicone layer may be aspirated and cultured.

If the silicone layer separates from the wound bed after neodermis formation begins, only the loose area of the silicone layer need be removed.

Either the wound should be grafted with an epidermal autograft immediately, or protected to prevent formation of granulation tissue.

Inflammation

INTEGRA Dermal Regeneration Template does not become inflamed unless there is a bacterial complication. This should be treated based on the clinical judgment of the practitioner.

HOW SUPPLIED

INTEGRA Dermal Regeneration Template is available in the following sizes:

- 5 cm x 5 cm (2 inch x 2 inch)
- 10 cm x 12.5 cm (4 inch x 5 inch)
- 10 cm x 25 cm (4 inch x 10 inch)
- 20 cm x 25 cm (8 inch x 10 inch)

The sheets consist of a porous matrix of cross-linked collagen and a glycosaminoglycan with an outer removable silicone covering identified by black sutures as markers to ensure proper placement on the wound bed.

The sheets of INTEGRA Dermal Regeneration Template Single Layer consist only of the dermal regeneration layer.

Each sheet of INTEGRA template is packaged in phosphate buffer in a foil pouch. Each foil pouch is packaged in a sealed outer chevron-style pouch.

Store flat at +2° to +30°C. Protect from freezing.

Once the product packaging has been opened, the product cannot be stored for reuse.

Please refer to the clinical training materials for complete instructions.

For product information or questions pertaining to sales and services, please contact the national distributor in your area or the manufacturer.

PRODUCT INFORMATION DISCLOSURE

INTEGRA LIFESCIENCES HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND THE MANUFACTURE OF THESE PRODUCTS. INTEGRA LIFESCIENCES EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA LIFESCIENCES SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA LIFESCIENCES NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS

SYMBOLS USED ON LABELING

	Do not re-use
	This product does not contain and is not manufactured with Dry Natural Rubber or Natural Rubber Latex
	Do not re-sterilize
	Do not use if package is damaged
	Expiration date
	Sterilized using irradiation
	Product complies with requirements of directive 93/42/EEC
	Lot number
	Consult Instructions for Use
	Temperature Limitation +2°C to +30°C
	Authorized Representative in the European Community
	Manufacturer
	Catalog Number
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner

INTEGRA™

 Integra LifeSciences Corporation
311 Enterprise Drive
Plainsboro, NJ 08536 USA

 Integra LifeSciences Services (France) SAS
Immeuble Sequoia 2
97 allée Alexandre Borodine
Parc Technologique de la Porte des Alpes
69800 Saint Priest – France
Telephone: 33 (0) 4 37 47 59 10
Fax: 33 (0) 4 37 47 59 29

Integra and Integra Dermal Regeneration Template are registered trademarks of Integra LifeSciences Corporation in the United States and/or other countries. Integra and the Integra logo are trademarks of Integra LifeSciences Corporation. ©2011 Integra LifeSciences Corporation. All rights reserved.

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